

SafeBasics® True Fit Thin

Nitrile Medical Examination Gloves

**Versatile,
all-purpose glove.**



■ **1185 A-E** (XS-XL)
 3 mils¹, 3.5 g
 300 gloves/box; 8 boxes/case

■ **1186 A-D** (XS-L)
 2.5 mils¹, 3.2 g
 300 gloves/box; 8 boxes/case

5+3



SafeTouch® Advanced Form

Nitrile Medical Examination Gloves

Protection and tactile sensitivity.

■ **1177 A-E** (XS-XL)
 2.8 mils¹, 3.7 g
 200 gloves/box; 10 boxes/case

8+2



AssureTouch® Basics

Nitrile Medical Examination Gloves

Dexterity and protection.

■ **9985 A-E** (XS-XL)
 2.7 mils¹, 3.5 g
 100 gloves/box; 10 boxes/case

7+3



SafeTouch® Connect

Latex Medical Examination Gloves

Comfort and grip.

■ **1124 A-E** (XS-XL)
 4.7 mils¹, 6.1 g
 100 gloves/box; 10 boxes/case

8+2



Ocean Pacific® NeoNatural

Polychloroprene Medical Examination Gloves

**Latex-free, strength of nitrile
with the comfort of latex.**

■ **1228 A-E** (XS-XL)
 4 mils¹, 6 g
 100 gloves/box; 10 boxes/case

7+3



AssureTouch® Guard Micro


Sterile Polychloroprene Surgical Gloves

The material shift to latex free.

Polychloroprene combines nitrile-like strength with latex-like comfort, without the risks of Type I allergies or latex sensitization. Ideal for procedures requiring tactile sensitivity. Micro thickness allows for double gloving.

■ **1144** (5.5-9)
 4.33 mils¹, 9.8 g
 50 pairs/box; 4 boxes/case

3+1



AssureTouch® Sense

Sterile Polychloroprene
Surgical Gloves

**Nitrile-like strength
with latex-like comfort.
Featuring 12" cuffs.**



3+1

■ **1145** (5.5-9)
 5.91 mils¹, 13.5 g
 50 pairs/box; 4 boxes/case

AssureTouch® **NEW** Defend

Sterile Polychloroprene
Surgical Gloves

**Nitrile-like strength
with latex-like comfort.
Featuring 12" cuffs.**



3+1

■ **1146** (5.5-9)
 5.91 mils¹, 13.5 g
 50 pairs/box; 4 boxes/case

AssureTouch® React

Sterile Polyisoprene
Surgical Gloves

**Latex-like elasticity
and comfort without
Type I allergy risks.**



3+1

■ **1147** (5.5-9)
 6.69 mils¹, 11.9 g
 50 pairs/box; 4 boxes/case

SafeMask® SofSkin

Procedure Earloop Face Masks

Ultra-soft and skin-friendly.



SafeMask® SofSkin
 Designed for sensitive skin.

ASTM LEVEL 3

■ 2086 ■ 2089

■ 2087 ■ 2088

50 masks/box, 10 boxes/case
 Also available in ASTM level 1.

3+1

SafeMask® SofSkin Fit **NEW**

Integrated chin piece
 provides a personalized fit.

ASTM LEVEL 3

■ 219914

50 masks/box, 10 boxes/case

3+1

SafeMask® SofSkin ClearFit **NEW**

Includes chin piece and
 anti-fog strip.

ASTM LEVEL 3

■ 220014

50 masks/box, 10 boxes/case

3+1

SafeMask® FreeFlow

Procedure Earloop Face Masks

Integrated strip prevents eyewear from fogging.



ASTM LEVEL 3

■ 200514 ■ 200511

■ 200515 ■ 200516

50 masks/box; 10 boxes/case.
 Also available in ASTM levels 1 & 2.

3+2

SafeMask® Classics

Procedure Earloop Face Masks

Reliable protection at a great price.



6+4

ASTM LEVEL 3

■ 205314

■ 205315

■ 205316

■ 205317

50 masks/box; 10 boxes/case
 Also available in ASTM levels 1 & 2.

SafeMask® Premier Elite

Procedure Earloop Face Masks

**Premium comfort.
Exceptional value.**

ASTM LEVEL 3

2047 2044 2043

2042 2046

50 masks/box; 10 boxes/case
 Also available in ASTM levels 1 & 2.

2+1



SafeMask® Artisan

Procedure Earloop Face Masks

**Innovative filtration layer provides
ultra lightweight feel.**



2+2

ASTM LEVEL 3

221414 221492 221415

50 masks/box; 10 boxes/case
 221415 also available in ASTM level 2.



SafeMask® Curve **NEW**

Procedure Earloop Face Masks

**Contoured fit ensures breathability and more
complete coverage.**



2+1

ASTM LEVEL 3

220114 220115

25 masks/box; 8 boxes/case

FREE FROM MOST COMMON ALLERGENS.*

Duraflor® Ultra

5% Sodium Fluoride Varnish

Superior uptake for lasting protection.

Non-clumping, effective treatment of hypersensitivity and cavity prevention.

- Dries clear and smooth
- Available in patient-friendly flavors

- | | |
|--|---|
|  1031-ST30 Strawberry |  1031-CH30 Cherry |
|  1031-M30 Mint |  1031-BG30 Bubble gum |

0.4 ml x 30 units/box. Promotion also applicable on 200 units/box

4+1


*Does not contain gluten, peanuts, tree nuts, soy, sesame, corn, shellfish, egg or dairy products.

DentiCare™ Swivel **NEW**

Oscillating Disposable Prophylaxis Angle

Splatter control?
Put a twist on it.


- | |
|--|
|  10057 Soft |
|  10058 Firm |
|  10059 Brush |

100 units/box

2+1

DentiCare™ Pro-Polish

Prophylaxis Paste with Fluoride

Low splatter and easy to rinse
in flavours patients will enjoy.


- | | |
|--|-------------------|
|  10047-FBGUN | Fine, Bubble Gum |
|  10047-MCHUN | Medium, Cherry |
|  10047-MMUN | Medium, Mint |
|  10047-CMUN | Coarse, Mint |
|  10047-CRBUN | Coarse, Raspberry |

 2 g/unit dose; 200/box
 Promotion also applicable on jars.

4+1

NEW

AssureCheck™ 10-Hour Biological Monitoring

Test. Incubate. Validate.

Buy 1 incubator
and get 1 free
box of ST1010

AssureCheck™ Dry Block Incubator
ST1000
1 unit/box

AssureCheck™ 10HR
Biological Indicator
ST1010
25 vials/box

4+1


SafeBasics® Sterilization Pouches

Convenient self-sealing design
with built-in internal and external
indicators at an affordable price.

99000 2 ¾ x 9"	99015 2 ¼ x 4"
99005 3 ½ x 5 ¼"	99025 5 ¼ x 10"
99010 3 ½ x 9"	99030 7 ½ x 13"

200 pouches/box

99010-B 3 ½ x 9"
500 pouches/box; 6 boxes/case

3+1


SafeSeal™ Quattro Sterilization Pouches

Self-Sealing

Engineered for peace of mind.

88000-4 2 ¾ x 9"	88025-4 5 ¼ x 10"
88005-4 3 ½ x 5 ¼"	88030-4 7 ½ x 13"
88010-4 3 ½ x 9"	88035-4 10 x 14"
88015-4 2 ¼ x 4"	88040-4 12 x 17"
88020-4 5 ¼ x 6 ½"	88045-4 4 ¼ x 11"

200 pouches/box

88010-4B 3 ½ x 9"
500 pouches/box; 6 boxes/case

2+1


Hedy® Quadra Fold SMMS Wrap

Four layers for exceptional
strength and protection.

620WR-15MW	15" x 15"
620WR-20MW	20" x 20"
620WR-24MW	24" x 24"

125 units/box

2+1


NEW

AssureCheck™ 20-Minute Biological Monitoring

**Fast, reliable results and
assured compliance.**

*Buy 1 incubator
and get 2 free
boxes of ST1020*

**AssureCheck™ Smart Incubator
ST2000**
1 unit/box

**AssureCheck™ Rapid
20 Minute Biological Indicator
ST1020**
50 vials/box

4+1


AssureCheck™ Type 5 **NEW**

Chemical Integrator

**Clear visual confirmation
that all parameters – time,
temperature, and presence
of steam – have been met.**

ST0005
250 units/bag

2+1


AssureCheck™ Multivariable **NEW**

Chemical Integrator

**Confirms proper time
and temperature conditions
have been met during
steam sterilization.**

(Type 4)
ST0004
250 units/bag

2+1
